


ACT-24

Automatic System for Blood Group Serology
The Only Universal and Standard Technique



- ***sedimentation analysis***
- ***image processing and visualization***



- ***internet communication***
- ***connection to hospital information system***
- ***special robotic dispenser is available***

Open system: spin tube technique

Applications:

- antibody screening
- antibody identification
- cross match
- ABO/RhD determination
- phenotyping
- titration studies

Methods:

- direct agglutination tests
- antiglobulin test (LISS, PEG)
- enzyme test
- polybrene test
- polybrene-AGT test

Automatic:

- cell washing
- reagent dispensing
- reading results
- printing results

Sample processing:

20-45 patients / hour (up to method)

CYTOTECH Ltd.

H-1131 Budapest, Dolmány u. 18-20

tel/fax: (+36)1-359-9360

www.cytotech.hu - info@cytotech.hu

 **CYTOTECH**
INSTRUMENTS LTD.